



IntelliPAP® AutoBilevel

Introducing the latest addition to the IntelliPAP line - the IntelliPAP AutoBilevel system incorporates the comfort of AutoAdjust® and Bilevel Technology to deliver comfortable therapy particularly at higher prescription pressures. The IntelliPAP AutoBilevel automatically adjusts to the user's pressure requirements on both inhalation and exhalation, using the unique DeVilbiss AutoAdjust algorithm featuring event set measurement.

Features Patented SmartCode for remote collection of Efficacy and Usage Data. Patent US 8649510



Pressure Settings	3-25 cm H ₂ O
While Breathing Compliance	Yes
Event Detection	Yes
Flow Rounding® Pressure Transition Relief	Yes
IPAP Rounding	Range: 0-5, increments of 1
EPAP Rounding	Range: 0-5, increments of 1
Fixed Pressure Support	Adjustable from 0-12 cm H ₂ O
Trigger Sensitivity	
Inspiratory	Range: 1-10, increments of 1
Expiratory	Range: 1-10, increments of 1
Automatic Leak Compensation	Yes
Remote Control Capabilities	Yes
Onboard Filter Clean Reminder	Yes
Onboard Compliance Quick Code	SmartCode®
Therapy Management Software	SmartLink®
Detailed Usage Data	SmartLink
Data Transfer via Memory Card	SmartLink
Remote Compliance Data Retrieval	SmartCode
Remote Efficacy Data Retrieval	SmartCode
Delay Pressure	Range: 3 cm H ₂ O to EPAP (or CPAP), increments of 0.5
Delay Time	0 to 45 minutes in 5-minute increments
Auto ON and Auto OFF	Yes
Visual Mask Off Alert	Yes
Product Weight	4.45 lbs with heated humidification system
Product Width	6.5 in.
Product Length/Depth	8.4 in.
Product Height	6.4 in.
Warranty (North America only)	3-year warranty on IntelliPAP 3-year warranty on heated humidification system 3-month warranty on humidifier chamber