

## FEELLIIFE Mesh Nebulizers



### Product description:

- 1.Portable design, small size, for home & office use.
- 2.Dimension:36mm\*36mm\*109mm, Weight:70g (without batteries).
- 3.The device adopting ultra-low noise design enables user feeling comfortable,so quite
- 4.Low voltage Indicator: the low voltage indicator light will turn red if battery weak, then the machine will stop work and power off automatically one minute later.
- 5.One-button operation design,high-output fine particle, Available to operate in virtually any position
- 6.Automatic shut-off after each medication dose (10 minutes use).
- 7.Charging via micro USB cable to any 5v usb power adapter or plug it directly into your computer.

### Product parameter:

Name	Portable Ultrasonic Nebulizer/Mesh Nebulizer
Weight	Approx. 70g (excluding batteries)
Medication Bottle Capacity	10ml maximum
Nebulization Rate	0.15ml/min~0.90ml/min
Power Consumption	<2W
Ultrasonic Frequency	110KHz±10KHz
Particle Size	1-5um

Sound Level	< 29dB
Power Supply	DC 7.4 V ( two AA lithium batteries) or 5.0V micro USB
Dimension	36mm*36mm*109mm

