

Oxygen Concentrators

9F-5AW Oxygen Concentrator

The 9F-5AW Oxygen Concentrator has a built-in pulse-oximeter with retractable cable and a large back-lit LCD screen for continuous monitoring of oxygen saturation of the patient. The back light has adjustable brightness for comfort. It has a nebulizer outlet, built-in timer and includes a remote control for easy control by the patient. The flow meter is also back-lit. This model features a self-diagnosis system and comes with multiple safety alarms. The precision compressor is made with German technology and has a service life of over 30,000 hours. The air inlet is situated at the top of the unit to reduce intake of dust particles and improve longevity.



Features:

- Built-in Pulse Oximeter
- Back-lit LCD screen
- Back-lit flow meter
- Remote control
- Nebulizer air outlet
- Top tray for storage of accessories
- Power failure alarm
- Low oxygen alarm
- Pressure alarm
- Compressor failure alarm

Technical Specifications:

Power : AC 220V @ 50 Hz, 400 VA

Flow range : 0.5 – 5.0 LPM

Oxygen concentration : 93% +/- 3%

Output pressure : 40-70 kPa

Operating noise : < 49.5 dB (front), 52 dB (overall)

Net weight : 18kg

